



ACCESSORIES

LIFEPAK® 15 MONITOR/DEFIBRILLATOR



GENUINE ACCESSORIES FROM PHYSIO-CONTROL

POWER OPTIONS

for the LIFEPAK 15 MONITOR/DEFIBRILLATOR

Thank you for choosing Physio-Control as your partner in helping you save lives and improve patient care. You already know our reputation for quality and reliability in LIFEPAK defibrillators and AEDs. When it comes to LIFEPAK accessories, we know you won't accept anything less. Inside this catalog you'll find accessories designed to meet your needs and reliability standards. All the accessories and disposables in this catalog are designed to work with the LIFEPAK 15 monitor/defibrillator. If you have any questions about accessories or disposables, please contact your Physio-Control representative.



Station Battery Charger

Includes AC power cord, mounting bracket and operating instructions. For use with the LIFEPAK 15 monitor/defibrillator. Batteries not included.

11577-000004

Mobile Battery Charger

Includes AC and DC power cords, mounting bracket and operating instructions. For use with the LIFEPAK 15 monitor/defibrillator. Batteries not included.

11577-000011



Lithium-ion Battery 5.7 Ah

5.7 amp hour, 11.1 volt, rechargeable, with fuel gauge. For use with the LIFEPAK 15 monitor/defibrillator.

21330-001176

ECG MONITORING ACCESSORIES

for the LIFEPAK 15 MONITOR/DEFIBRILLATOR



5-Wire ECG Cable
Right-angle connector, 4-wire limb plus 1 chest lead, labeled "V1" on the LIFEPAK 15 monitor reports.
11110-000066



3-Wire ECG Cable
Right-angle connector
11110-000029



4-Wire Limb Lead Attachment
11110-000102



12-Lead ECG Patient Cable
Main 5-foot trunk, right-angle connector
11110-000111

Main 8-Foot Trunk, right-angle connector
11110-000110



6-Wire Precordial Lead Attachment
11110-000103



LIFE•PATCH® ECG Electrodes
Adult, pregelled (package of 3 electrodes)
11100-000001



Box of Strip Chart Recorder Paper
100mm x 22m (2 rolls per box)
11240-000016



LIFE•PATCH ECG Electrodes
Adult, pregelled (package of 4 electrodes)
11100-000002

THERAPY DELIVERY ACCESSORIES

for the LIFEPAK 15 MONITOR/DEFIBRILLATOR

HARD PADDLES AND ELECTRODE GEL



Standard Hard Paddles English

1 pair, for use with the LIFEPAK 15 monitor/defibrillator
11130-000061



Pediatric Paddle, External

One each (two required); Slips onto standard adult hard paddle.
11133-000001



SIGNAGEL® Electrode Gel

One (1) 8.5 oz. tube for use with hard paddles. Highly conductive, multi-purpose electrolyte meets all the standards of the ideal saline electrode gel. Recommended for ECG, defibrillation, biofeedback and EMG.

21300-005847

EDGE SYSTEM™ ELECTRODES FOR PACING/DEFIBRILLATION/ECG WITH QUIK-COMBO® CONNECTOR

18-month minimum shelf life remaining at time of shipment from Physio-Control except where noted.



EDGE System Electrodes with QUIK-COMBO Connector

24" leadwire length
11996-000091



EDGE System RTS (Radiotransparent) Electrodes with QUIK-COMBO Connector

24" leadwire length
11996-000090



EDGE System Electrodes with QUIK-COMBO Connector and REDI-PAK™ Preconnect System

42" leadwire length
11996-000017



Pediatric EDGE System RTS Electrodes with QUIK-COMBO Connector

For use only with manual monitor/defibrillators; 12 month minimum shelf life at time of shipment
24" leadwire length
11996-000093



QUIK-COMBO® Therapy Cable

With convenient TRUE-LOCK™ Cable Connector. Length is approximately 8 ft
11113-000004

SpO₂, SpCO AND SpMet MONITORING ACCESSORIES

for the LIFEPAK 15 MONITOR/DEFIBRILLATOR

MASIMO® SET® RAINBOW® DIRECT CONNECT SENSORS



Masimo SET Rainbow DCI-dc3, Adult Reusable Direct Connect Sensor

3 ft Reusable direct connect cable and sensor for patients > 30 kg, for use with LIFEPAK 15 monitor/defibrillator

11996-000335

Masimo SET Rainbow DCI-dc8, Adult Reusable Direct Connect Sensor,

8 ft Reusable direct connect cable and sensor for patients >30 kg, for use with LIFEPAK 15 monitor/defibrillator

11171-000032

Masimo SET Rainbow DCI-dc12, Adult Reusable Direct Connect Sensor

12 ft Reusable direct connect cable and sensor for patients > 30 kg, for use with LIFEPAK 15 monitor/defibrillator

11996-000336

Masimo SET Rainbow DCIP-dc3, Pediatric Reusable Direct Connect Sensor

3 ft reusable direct connect cable and sensor for patients 10-50 kg, for use with LIFEPAK 15 monitor/defibrillator

11996-000337

Masimo SET Rainbow DCIP-dc8, Pediatric Reusable Direct Connect Sensor

8 ft reusable direct connect cable and sensor for patients 10-50 kg, for use with LIFEPAK 15 monitor/defibrillator

11171-000033

Masimo SET Rainbow DCIP-dc12, Pediatric Reusable Direct Connect Sensor

12 ft reusable direct connect cable and sensor for patients 10-50 kg, for use with LIFEPAK 15 monitor/defibrillator

11996-000338

SpO₂ MONITORING ACCESSORIES

for the LIFEPAK 15 MONITOR/DEFIBRILLATOR

MASIMO SET LNCS[®] SENSORS AND CABLES



Masimo SET Red LNC-04 LNCS Patient Cable

4 ft reusable connector cable, for use with LIFEPAK 15 monitor/defibrillator
11996-000323

Masimo SET Red LNC-10 LNCS Patient Cable

10 ft reusable connector cable, for use with LIFEPAK 15 monitor/defibrillator
11996-000324

Masimo SET Red LNC-14 LNCS Patient Cable

14 ft reusable connector cable, for use with LIFEPAK 15 monitor/defibrillator
11996-000325



Masimo SET LNCS DCI Adult Reusable Sensor

3 ft multiuse sensor for patients >30 kg, for use with LNCS patient cables on LIFEPAK 12 or LIFEPAK 15 monitor/defibrillators.
11171-000017



Masimo SET LNCS Pdtx Pediatric Adhesive Sensor

18" single patient use sensor for patients 10–50 kg (20/box), for use with LNCS patient cables on LIFEPAK 12 or LIFEPAK 15 monitor/defibrillators
11171-000020



Masimo SET LNCS Amtx Adult Adhesive Sensor

18" single patient use sensor for patients >30 kg (20/box), for use with LNCS patient cables on LIFEPAK 12 or LIFEPAK 15 monitor/defibrillators
11171-000019



Masimo SET LNCS DCIP Pediatric Reusable Sensor

3 ft multiuse sensor for patients 10–50 kg (box of 20), for use with LNCS patient cables on LIFEPAK 12 or LIFEPAK 15 monitor/defibrillators
11171-000018

SpO₂ MONITORING ACCESSORIES

for the LIFEPAK 15 MONITOR/DEFIBRILLATOR

MASIMO SET LNCS SENSORS AND CABLES (CONTINUED)



Masimo SET LNCS Neo-L Neonatal Adhesive Sensor

18" single patient use sensor for patients <3 kg or >40 kg (20/box), for use with LNCS patient cables on LIFEPAK 12 or LIFEPAK 15 monitor/defibrillators

11171-000028



Masimo SET LNCS NeoPt-L Neonatal Adhesive Sensor

18" single patient use sensor for patients <1 kg (20/box), for use with LNCS patient cables on LIFEPAK 12 or LIFEPAK 15 monitor/defibrillators.

11171-000029



Masimo SET LNCS Inf-L Infant Adhesive Sensor

18" single patient use sensor for patients 3–20 kg (20/box), for use with LNCS patient cables on LIFEPAK 12 or LIFEPAK 15 monitor/defibrillators

11171-000031

MASIMO SET RED™ DIRECT CONNECT SENSORS



Masimo SET Red DCI-dc3 Adult Reusable Direct Connect Sensor

3 ft Reusable direct connect cable and sensor for patients >30 kg, for use with LIFEPAK 15 monitor/defibrillator

11996-000331

Masimo SET Red DCI-dc12 Adult Reusable Direct Connect Sensor

12 ft Reusable direct connect cable and sensor for patients >30 kg, for use with LIFEPAK 15 monitor/defibrillator

11996-000332



Masimo SET Red DCI-dc3 Pediatric Reusable Direct Connect Sensor

3 ft Reusable direct connect cable and sensor for patients 10-50 kg, for use with lifePaK 15 monitor/defibrillator

11996-000333

Masimo SET Red DCIP-dc12, Pediatric Reusable Direct Connect Sensor

12 ft Reusable direct connect cable and sensor for patients 10-50 kg, for use with LIFEPAK 15 monitor/defibrillator

11996-000334

END-TIDAL CO₂ (ETCO₂) MONITORING ACCESSORIES

for the LIFEPAK 15 MONITOR/DEFIBRILLATOR

ORIDION® INTUBATED FILTERLINES



FilterLine® SET Adult/Pediatric
Includes airway adapter. FilterLine (79") for short-term intubated patients, 24 hours typical.
11996-000299 (pack of 10)
11996-000081 (box of 25)



FilterLine H SET Adult/Pediatric
Includes airway adapter, High-humidity FilterLine (79") for long-term intubated patients, 72 hours typical (box of 25).
11996-000080

FilterLine SET Long Adult/Pediatric

Includes airway adapter, Long-length FilterLine (185") for short-term intubated patients, 24 hours typical (box of 25).
11996-000164



FilterLine H SET Infant/Neonatal
Includes airway adapter, High-humidity FilterLine (79") for long-term intubated patients, 72 hours typical (box of 25).
11996-000001

ORIDION NON-INTUBATED FILTERLINES



Smart CapnoLine® Plus with O₂ Delivery Adult/Intermediate
For patients above 44 lbs, O₂ delivery and CO₂ sampling via mouth or nose, 12 hours equals typical length of use.
11996-000298 (pack of 10)
11996-000163 (box of 25)
11996-000167 (case of 100)



Smart CapnoLine with O₂ Delivery Pediatric
For patients less than 44 lbs, O₂ delivery and CO₂ sampling via mouth or nose, 12 hours equals typical length of use (box of 25).
11996-000128

Smart CapnoLine Plus Long with O₂ Delivery Adult/Intermediate

For patients above 44 lbs, long-length nasal cannula (165") for O₂ delivery and CO₂ sampling via mouth or nose, 12 hours equals typical length of use (box of 25).
11996-000165



Smart CapnoLine Pediatric
For patients less than 44 lbs, CO₂ sampling via the mouth or the nose, 12 hours equals typical length of use (box of 25).
11996-000120

Smart CapnoLine Plus Adult/Intermediate

For patients above 44 lbs, CO₂ sampling via the mouth or the nose, 12 hours equals typical length of use.
11996-000297 (pack of 10)
11996-000162 (box of 25)
11996-000166 (case of 100)

NIBP MONITORING ACCESSORIES

for the LIFEPAK 15 MONITOR/DEFIBRILLATOR



Adult Cuff
14 x 37.5 cm, reusable
11996-000024

Child Cuff (Pediatric)
9 x 27 cm, reusable
11996-000022

Small Adult Cuff
12 x 30 cm, disposable
11996-000029

Child Cuff (Pediatric)
9 x 22.5 cm, disposable
11996-000028

Large Adult Cuff
16 x 42 cm, reusable
11996-000025

Small Child Cuff
7 x 21 cm, reusable
11996-000021

Adult Cuff
14 x 35 cm, disposable
11996-000030

Infant Cuff
6 x 15 cm, disposable
11996-000027

Small Adult Cuff
12 x 30 cm, reusable
11996-000023

Thigh Cuff
20 x 54 cm, reusable
11996-000026

Large Adult Cuff
16 x 40 cm, disposable
11996-000031

Infant Cuff
6 x 18 cm, reusable
11996-000020

X-Large Adult Cuff
20 x 50 cm, disposable
11996-000032



NIBP Tubing
Coiled, 2-9 feet, for use with
LIFEPAK 15 monitor/defibrillator
21300-007300



NIBP Tubing
9 feet, for use with
LIFEPAK 15 monitor/defibrillator
21300-007299

12 feet, for use with
LIFEPAK 15 monitor/defibrillator
21300-007298

CARRYING CASES, COMMUNICATION ACCESSORIES AND MATERIALS

for the LIFEPAK 15 MONITOR/DEFIBRILLATOR



Standard Carrying Case

Includes right pouch and left pouch, for use with LIFEPAK 15 monitor/defibrillator
11577-00002

Shoulder Strap

For use with LIFEPAK 15 monitor/defibrillator
11577-00001



Back Pouch

Ideal for additional accessory storage, for use with LIFEPAK 15 monitor/defibrillator
11260-000039

Top Pouch

Storage for sensors and electrodes; insert in place of standard paddles
11220-000028



LIFEPAK Monitor to Serial Port Cable

For connecting LIFEPAK 12 or LIFEPAK 15 monitor/defibrillator to PC
11230-000020



LIFEPAK Monitor to PC USB Cable

For connecting LIFEPAK 12 or LIFEPAK 15 monitor/defibrillator to PC
11996-000288



Analog ECG Output Cable

11110-000044

TRAINING TOOLS AND TESTERS

for the LIFEPAK 15 MONITOR/DEFIBRILLATOR



QUIK-COMBO 12-Lead Patient Simulator

Connects directly to your LIFEPAK defibrillator for safe simulation of cardioversion and electrical capture. Generates fibrillation, tachycardias, and bradycardias, as well as ST segment and T wave abnormalities. For use with LIFEPAK devices with a 12-lead ECG feature.

11996-000311

QUIK-COMBO 3-Lead Patient Simulator

Connects directly to your LIFEPAK defibrillator for safe, interactive training. Select from 17 ECG rhythms including: fibrillation, tachycardias and bradycardias.

11996-000310



Clip-on Training Electrodes

For use with QUIK-COMBO patient simulators. Allows users to practice defibrillation electrode placement. Includes instructions for use and five (5) pair of reusable training electrodes.

11250-000052

Replacement Training Electrodes

Includes five (5) pair of reusable training electrodes. Adhesive electrode portion only.

11101-000003

Operating Instructions: LIFEPAK 15 Monitor/Defibrillator

Also available as free online download at www.physio-control.com

26500-002408

Operating Instructions: Station and Mobile Battery Charger

26500-002948

Inservice Video: LIFEPAK 15 Monitor/Defibrillator

DVD format in English

21330-001357



12-Leads Made Easy

Web-based training program covering 12-Lead ECG basics, STEMI recognition, STEMI Imposters, rapid 12-lead acquisition, and over 125 real 12-lead ECG examples. Eligible EMS and nursing participants receive 5 hours of continuing education credit upon successful completion of the final exam. Price is per person for one year.

44500-000001



Defibrillator Checker

Tests integrity of energy delivery through standard hard paddles. Neon light indicates energy has been delivered.

11998-000060

Test Load

Use to perform therapy cable performance checks. Connects to the QUIK-COMBO therapy cable on the defibrillator

21330-001365 (English only)

21330-001367 (English and French)



To order accessories for the LIFEPAK 15 monitor/defibrillator, visit our webstore at <https://store.physio-control.com> (US only) or contact us at one of the numbers listed below



Physio-Control Headquarters
11811 Willows Road NE
Redmond, WA 98052
www.physio-control.com

Customer Support
P. O. Box 97006
Redmond, WA 98073
Toll Free 800 442 1142
Fax 800 426 8049

Physio-Control Canada
Medtronic of Canada Ltd
99 Hereford Street
Brampton, ON
L6Y 0R3
Tel 888 879 0977
Fax 866 430 6115

