



## MIDMARK 630 PROCEDURE TABLE

# FEATURES

No matter what type of procedures you perform, proper positioning is crucial to providing accessibility and flexibility at the point of care. Power height, back, foot, and tilt gives the Midmark 630 virtually endless movement capabilities, allowing you to position patients correctly for the most effective procedure. Adjustable height allows you to work in a position that is comfortable for you.

Procedures can be stressful and intimidating to patients. The Midmark 630 is designed to position patients for maximum comfort. Added features like the patented Premium Comfort System and optional heated upholstery helps ensure a calmer patient, reducing movement during procedures.



**CeviMed**  
MEDICAL EQUIPMENT & SUPPLIES

## MIDMARK 630 PROCEDURE TABLE

### SPECIFICATIONS

Patient Weight (Maximum).....450 lbs  
(204.1 kgs)  
Paper Roll (maximum size): .....18 in. long  
x 3.5 in. diameter  
(45.7 x 8.9 cm)  
Range of Motion &  
Dimensions .....[See illustration]  
Weight of Table:  
230 / 630 / 75 (w/standard base) .....455 lbs  
(206.4 kg)  
w/packaging & skid .....510 lbs  
(231.3 kg)  
230 / 630 / 75 (w/rotational base) .....540 lbs  
(244.9 kg)  
w/packaging & skid .....600 lbs  
(271.1 kg)  
Power Cord: .....8 ft. (243.8  
cm) long  
Electrical Requirements: .....[See  
Compliance Charts]  
Foot / Hand Control Voltage: .....10 VAC,  
SELV (Safety Extra Low Voltage)  
Duplex Receptacle (optional)  
maximum load: .....115 VAC,  
3A, 50/60 Hz  
Fuses:  
Main Power (located at power cord inlet):  
F1 & F2 .....6.3A, 250V,  
Type-T, 5 x 20 mm



1 year warranty.



Brand new upholstery.



Save 20% - 50% on  
Procedure tables by  
buying refurbished.

